
MediWet Humidifier. Product description.



Product description

The humidifier for oxygen therapy is a device allowing an increase in the relative humidity in the oxygen supplied to the patient; this gas for medical use normally lacks a sufficient humidity grade to be physiologically tolerated.

The jar is made of polycarbonate, the cap with inlet adaptor of polypropylene and inlet connection of polypropylene + coated brass. The oxygen therapy humidifier is autoclavable at 121°C (15 min).

Technical data

Inner volume	200 ml
Height	190 ml
Diameter	57 mm
Width with inlet adaptor	67 mm
Weight	70 g

Article numbers

Art. No.	Inlet Connection	Outlet Connection
330385	9/16 UNF	Ø 5 - 9
330386	G3/8	Ø 5 - 9
330387	Quick connector MC 315	Ø 5 - 9

Spare parts

Art. No.	Description
330976	Silicone hose 10 pcs
330973	Silicone hose + diffuser 10 pcs
330484	Diffuser 10 pcs
330970	Washers 10 pcs (under lid)
331970	O-ring G3/8 inlet 10 pcs
331971	O-ring 9/16 UNF inlet 10 pcs

Recycling

The product shall be recycled according to local regulations.

Contact us Linde Healthcare is committed to working with healthcare providers and regulatory authorities to continuously promote safe use of medical products and improve patient care.

We provide medical gas products, therapies, technical solutions and services to hospitals, clinics, nursing facilities, emergency management services and home healthcare providers around the world.

With our long experience and understanding of healthcare realities, you can depend on solutions that are delivered and serviced to the highest possible standards of quality, safety and efficiency.